

Antibiotic Hemostasis Pad

ez CLOT Radial

Sterilized pad fabricated from 100% Chitosan fiber to stop bleeding from the radial artery after the treatment.



How to use it.

1. After pulling the sheath out by about 3 -5 cm, put ezCLOT R on the puncture site.
2. Fix ezCLOT R with a surgical bandage.
3. While pressing ezCLOT R with a thumb, remove the sheath.
4. Bind ezCLOT R with a cohesive bandage at most 4 or 5 times for pressure.
5. If a patient suffers from the pressure, unwind the cohesive bandage 1 or 2 times to the opposite direction and cut it unwound.
6. Remove it after confirming bleeding stop within 2 hours.

Femoral
ez CLOT

Specification

	2.5 X 2.5 X 0.4 cm	1 Box : 10개입
	3 X 3 X 0.4 cm	1 Box : 10개입
	5 X 5 X 0.4 cm	1 Box : 10개입

Headquarter

B-1409, Woolim Lionsvalley 5th, 302 Galmachi ro, Joongwon gu, Seongnam si, Gyeonggi do, Korea / TEL +82 31 755 7302 / FAX +82 31 755 7306

Factory

1-207, Wonju Venture Factory, 130-2 Donghwagongdan ro, Munmak eup, Wonju si, Gangwon do, Korea / TEL +82 33 743 7302 / FAX +82 33 743 7304